

CERTIFICATE OF ANALYSIS

SERVA Triple Color Protein Standard III	Cat.No. : 39258
5 – 245 kDa, Liquid Mix	Lot No. : 130929

Protein Components	Molecular Weight (kDa)
Protein 1 (blue)	245
Protein 2 (blue)	180
Protein 3 (blue)	135
Protein 4 (blue)	100
Protein 5 (orange)	75
Protein 6 (blue)	63
Protein 7 (blue)	48
Protein 8 (blue)	35
Protein 9 (green)	25
Protein 10 (blue)	20
Protein 11 (blue)	17
Protein 12 (blue)	11
Protein 13 (blue)	5
Storage (°C)	-15 to -25
Minimum Shelf Life	08/ 2015

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.08.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

This report has been computer generated and does not contain a signature.