

CERTIFICATE OF ANALYSIS**SERVA Triple Color Protein Standard III**
5 – 245 kDa, Liquid Mix**Cat.No. : 39258**
Lot No. : 160950

Protein Components	Molecular Weight (kDa)
Protein 1 (blue)	245
Protein 2 (blue)	180
Protein 3 (blue)	135
Protein 4 (blue)	100
Protein 5 (orange)	75
Protein 6 (blue)	63
Protein 7 (blue)	48
Protein 8 (blue)	35
Protein 9 (green)	25
Protein 10 (blue)	20
Protein 11 (blue)	17
Protein 12 (blue)	11
Protein 13 (blue)	5
Storage (°C)	-15 to -25
Minimum Shelf Life	06.10.2018

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control**Printing date: 17.10.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.