

CERTIFICATE OF ANALYSIS

SERVA Triple Color Protein Standard III
5 – 245 kDa, Liquid Mix

Cat.No. : 39258
Lot No. : 230591

| Protein Components | Molecular Weight (kDa) |
|---------------------------|-----------------------------------|
| Protein 1 (blue) | 245 |
| Protein 2 (blue) | 180 |
| Protein 3 (blue) | 135 |
| Protein 4 (blue) | 100 |
| Protein 5 (orange) | 75 |
| Protein 6 (blue) | 63 |
| Protein 7 (blue) | 48 |
| Protein 8 (blue) | 35 |
| Protein 9 (green) | 25 |
| Protein 10 (blue) | 20 |
| Protein 11 (blue) | 17 |
| Protein 12 (blue) | 11 |
| Protein 13 (blue) | 5 |
| Storage (°C) | -15 to -25 |
| Minimum Shelf Life | 07.11.2025 |

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 05.06.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.