

**CERTIFICATE OF ANALYSIS**

<b>SERVA VisiBlot Standard I</b>	<b>Cat.No. : 39260</b>
25 – 150 kDa, Liquid Mix	<b>Lot No. : 180184</b>

<b>Protein Components</b>	<b>Molecular Weight (kDa)</b>
<b>Protein 1</b>	150
<b>Protein 2</b>	120
<b>Protein 3</b>	90
<b>Protein 4 (orange)</b>	85
<b>Protein 5</b>	75
<b>Protein 6</b>	60
<b>Protein 7 (orange)</b>	45
<b>Protein 8</b>	40
<b>Protein 9 (blue)</b>	25
<b>Protein 10</b>	25
<b>Storage (°C)</b>	-15 to -25
<b>Minimum Shelf Life</b>	02.05.2020

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 21.02.2018**

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.