

Certificate of Analysis

SERVA Fluo-610 Standard I**Cat. no. : 39261****Lot No.: P210270**

Protein Components	Molecular weight
Phosphorylase B	97 400
Albumin bovine (BSA)	67 000
Albumin egg	45 000
Carbonic anhydrase	29 000
Trypsin inhibitor (soybean)	21 000
Lysozyme (chicken egg white)	14 400
Minimum shelf life	26.01.2023
Storage (°C)	-15 to -25

Reconstitution

Dissolve content of one vial in 50 µl dist. water, mix carefully and spin down.

A volume of 5 µl is sufficient for SDS-PAGE (in a standard form of 80 x 60 x 0.75 mm).

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control**Printing date: 29.07.2021**

Dr. Judith Koch

Patricia Milford

This report has been computer generated and does not contain a signature.