

# Certificate of Analysis

---

**SERVA Fluo-610 Standard I****Cat. no. : 39261****Lot No.: P230172**

---

<b>Protein Components</b>	<b>Molecular weight</b>
<b>Phosphorylase B</b>	97 400
<b>Albumin bovine (BSA)</b>	67 000
<b>Albumin egg</b>	45 000
<b>Carbonic anhydrase</b>	29 000
<b>Trypsin inhibitor (soybean)</b>	21 000
<b>Lysozyme (chicken egg white)</b>	14 400
<b>Minimum shelf life</b>	27.04.2024
<b>Storage (°C)</b>	-15 to -25

**Reconstitution**

Dissolve content of one vial in 50 µl dist. water, mix carefully and spin down.

A volume of 5 µl is sufficient for SDS-PAGE (in a standard form of 80 x 60 x 0.75 mm).

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

---

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date: 02.05.2023**

Dr. Judith Koch

Marina Konter

This report has been computer generated and does not contain a signature.