

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

SERVA FastLoad 50 bp DNA Leiter	Kat.Nr./Cat.No. :	39315
SERVA FastLoad 50 bp DNA ladder	Lot/Contr.No. :	151391

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Gebrauchsfertiger DNA-Standard für Agarosegele. Besteht aus 17 Fragmenten im Bereich 50 Bp bis 1500 Bp. Ready-to-Use DNA standard for agarose gel electrophoresis containing 17 fragments from 50 – 1500 base pairs.		
Anwendungstest application test	Elektrophorese Electrophoresis	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			09.11.2020
Lagerung (°C) storage			-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.05.2019

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.