

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>SERVA FastLoad 50 bp DNA Leiter</b>	<b>Kat.Nr./Cat.No. :</b>	<b>39315</b>
SERVA FastLoad 50 bp DNA ladder	<b>Lot/Contr.No. :</b>	<b>240965</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Beschreibung</b> description	Gebrauchsfertiger DNA-Standard für Agarosegele. Besteht aus 17 Fragmenten im Bereich 50 Bp bis 1500 Bp. Ready-to-Use DNA standard for agarose gel electrophoresis containing 17 fragments from 50 – 1500 base pairs.		
<b>Anwendungstest</b> application test	Elektrophorese Electrophoresis	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			28.10.2026
<b>Lagerung (°C)</b> storage			-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 11.11.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.