

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS



<b>SERVA FastLoad 100 bp DNA Leiter</b>	<b>Kat.Nr./Cat.No. :</b>	<b>39316</b>
SERVA FastLoad 100 bp DNA ladder	<b>Lot/Contr.No. :</b>	<b>240505</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Beschreibung</b> description	Gebrauchsfertiger DNA-Standard für Agarosegele. Besteht aus 12 Fragmenten im Bereich 100 Bp bis 3000 Bp. Ready-to-Use DNA standard for agarose gel electrophoresis containing 12 fragments from 100 – 3000 base pairs.		
<b>Anwendungstest</b> application test	Elektrophorese Electrophoresis	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			03.06.2026
<b>Lagerung (°C)</b> storage			-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 11.06.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.