

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Denhardt's Solution, 50x Konzentrat	Kat.Nr./Cat.No. : 39603
Denhardt's solution, 50x concentrate	Lot/Contr.No. : P130376
molecular biology grade	

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare farblose bis leicht gelbliche Lösung clear colorless to slightly yellowish solution	entspricht corresponds
pH		7.3
Zusammensetzung composition		
Albumin Fraktion V	0.1g/ 10 ml	entspricht corresponds
Polyvinylpyrrolidone	0.1g/ 10 ml	entspricht corresponds
Ficoll [®] 400	0.1g/ 10 ml	entspricht corresponds
DNase/ RNase Aktivität DNase/ RNase activity	nicht nachweisbar none detected	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		19.09.2017
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.07.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.