

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|---|--------------------------------|
| Denhardt´s Solution, 50x Konzentrat | Kat.Nr./Cat.No. : 39603 |
| Denhardt´s solution, 50x concentrate molecular biology grade | Lot/Contr.No. : P160392 |

| Parameter parameter | Spezifikation specification | Ergebnis result |
|--|--|---------------------------|
| Aussehen appearance | klare farblose bis leicht gelbliche Lösung clear colorless to slightly yellowish solution | entspricht corresponds |
| pH | | 7.4 |
| Zusammensetzung composition | | |
| Albumin Fraktion V | 0.1g/ 10 ml | entspricht corresponds |
| Polyvinylpyrrolidone | 0.1g/ 10 ml | entspricht corresponds |
| Ficoll [®] 400 | 0.1g/ 10 ml | entspricht corresponds |
| DNase/ RNase Aktivität DNase/ RNase activity | nicht nachweisbar none detected | entspricht corresponds |
| Mindesthaltbarkeit minimum shelf life | | 19.10.2018 |
| Lagerung (°C) storage | | -15 bis/ to -25 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.10.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.