

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Tris-Puffer pH 8.8, 1M Lösung</b>	<b>Kat.Nr./Cat.No. : 39794</b>
Tris buffer pH 8.8, solution 1M	<b>Lot/Contr.No. : P160258</b>
molecular biology grade	

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
<b>Zusammensetzung</b> composition			
TRIS (#37186)		121.14 g/l	entspricht corresponds
<b>pH</b>	20 °C eingestellt mit HCl adjusted with HCl	8.7 – 8.9	8.9
<b>DNase/ RNase Aktivität</b> DNase/ RNase activity		nicht nachweisbar none detected	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			24.06.2018
<b>Lagerung (°C)</b> storage			+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 28.06.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.