

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Dowex[®] 1X2 (200 – 400 mesh) , pract.	Kat.Nr./Cat.No.: 41030 Lot/Contr.No.: 160181
---	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis hellelfenbein- farbenes Korn white to light-ivory colored particles	entspricht corresponds
Trocknungsverlust (%) loss on drying	3h, 105°C	70.0 – 80.0	71.1
Gegenion counter ion		Cl ⁻	entspricht corresponds
Kapazität (eq/l) capacity		min. 0.6	0.9
Mindesthaltbarkeit minimum shelf life			17.02.2021
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.02.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.