

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Dowex® 1X2 p.A.</b> (200 – 400 mesh) analytical grade	<b>Kat.Nr./Cat.No. : 41031</b> <b>Lot/Contr.No. : P120123</b>
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		weißes bis hellelfenbein- farbenes Korn white to light-ivory colored particles	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	105 °C, 3 h	70.0 – 80.0	76.8
<b>Gegenion</b> counter ion		Cl <sup>-</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 0.6	0.6
<b>Mindesthaltbarkeit</b> minimum shelf life			03/ 2017
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 26.03.2012**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.