

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>Dowex® 1X2 p.A.</b> (200 – 400 mesh) analytical grade	<b>Kat.Nr./Cat.No. :</b>	<b>41031</b>
	<b>Lot/Contr.No. :</b>	<b>P160009</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		weißes bis hellelfenbein- farbenes Korn white to light-ivory colored particles	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	105 °C, 3 h	70.0 – 80.0	70.1
<b>Gegenion</b> counter ion		Cl <sup>-</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 0.6	0.9
<b>Mindesthaltbarkeit</b> minimum shelf life			13.01.2027
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**
**Druckdatum/ printing date: 16.05.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.