

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>DOWEX® 1X8.</b> (20 – 50 mesh), pract.	<b>Kat.Nr./Cat.No. :</b> <b>41080</b>
	<b>Lot/Contr.No. :</b> <b>220625</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		hellelfenbeinfarbenes bis beiges Korn light ivory-coloured to beige bead	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	105 °C, 3 h	45.0 – 60.0	52.8
<b>Gegenion</b> counter-ion		Cl <sup>-</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 1.3	1.41
<b>Mindesthaltbarkeit</b> minimum shelf life			22.06.2027
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Dieser Artikel kann produktionsbedingt andersfarbige Partikel enthalten.**  
**Diese haben keine Auswirkungen auf die Eigenschaften und die Funktionalität des Produktes.**  
 Due to the production process this article can contain differently coloured particles.  
 This does not affect the quality characteristics and the functionality of the product.

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 04.07.2022**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.