

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

DOWEX® 1X8. (20 – 50 mesh), pract.	Kat.Nr./Cat.No. : 41080
	Lot/Contr.No. : 240650

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellelfenbeinfarbenes bis beiges Korn light ivory-coloured to beige bead	entspricht corresponds
Trocknungsverlust (%) loss on drying	105 °C, 3 h	45.0 – 60.0	52.6
Gegenion counter-ion		Cl ⁻	entspricht corresponds
Kapazität (eq/l) capacity		min. 1.3	1.40
Mindesthaltbarkeit minimum shelf life			23.07.2029
Lagerung (°C) storage			+15 bis/ to +30

**Dieser Artikel kann produktionsbedingt andersfarbige Partikel enthalten.
 Diese haben keine Auswirkungen auf die Eigenschaften und die Funktionalität des Produktes.**
 Due to the production process this article can contain differently coloured particles.
 This does not affect the quality characteristics and the functionality of the product.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 05.08.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.