

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>DOWEX® 1X8</b> , p.A. (20 – 50 mesh) analytical grade	<b>Kat.Nr./Cat.No. :</b> 41081 <b>Lot/Contr.No. :</b> P190371
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		sandgelbes, transparentes Korn sandyellow-colored, transparent particles	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	3 h, 105 °C	50.0 – 60.0	52.5
<b>Gegenion</b> counter-ion		Cl <sup>-</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 1.3	1.38
<b>Mindesthaltbarkeit</b> minimum shelf life			07.11.2024
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Das Produkt kann schwarze Partikel enthalten. Dabei handelt es sich um inaktive Copolymere, die die spezifizierten Eigenschaften wie Kapazität, Wassergehalt oder die Abgabe von organischen Substanzen vom Harz nicht beeinträchtigen.**

The product can contain black particles. These are inactive copolymers which do not impair the specified characteristics as capacity, water content or the release of organic substances of the resin.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 08.11.2019

Christian Monsler

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.