

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Dowex<sup>®</sup> 1X8 p.A.</b> (50 – 100 mesh) analytical grade	<b>Kat.Nr./Cat.No.:</b> 41091
	<b>Lot/Contr.No.:</b> P070193

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		hellelfenbeinfarbenes Korn light-ivory colored particles	entspricht corresponds
<b>Gegenion</b> counter ion		Cl <sup>-</sup>	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	105°C, 3h	45.0 – 55.0	46.9
<b>Kapazität (eq/l)</b> capacity		min. 1.2	1.6
<b>Mindesthaltbarkeit</b> minimum shelf life			07/ 2016
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 02.06.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.