

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

DOWEX® 1X8 p.A. (100 – 200 mesh) analytical grade	Kat.Nr./Cat.No. : 41101 Lot/Contr.No. : P140564
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		elfenbeinfarbenes Korn ivory colored particles	entspricht corresponds
Gegenion counter ion		Cl ⁻	entspricht corresponds
Trocknungsverlust (%) loss on drying	3 h, 105 °C	39.0 – 50.0	44.0
Kapazität (eq/l) capacity		min. 1.2	1.4
Mindesthaltbarkeit minimum shelf life			12/ 2019
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 17.12.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.