

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

DOWEX® 1X8 p.A. (100 – 200 mesh) analytical grade	Kat.Nr./Cat.No. : 41101 Lot/Contr.No. : P210326
-------------------------------------------------------------	------------------------------------------------------------------

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		elfenbeinfarbenes Korn ivory colored particles	entspricht corresponds
Gegenion counter ion		Cl ⁻	entspricht corresponds
Trocknungsverlust (%) loss on drying	3 h, 105 °C	39.0 – 50.0	42.3
Kapazität (eq/l) capacity		min. 1.2	1.4
Mindesthaltbarkeit minimum shelf life			24.09.2026
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.10.2021

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.