

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>DOWEX® 50 WX2</b> p.A. (100 – 200 mesh), analytical grade	<b>Kat.Nr./Cat.No. :</b> 41521 <b>Lot/Contr.No. :</b> P130454
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		elfenbeinfarbenes Korn ivory-colored particles	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	3 h, 105 °C	74.0 – 82.0	80.5
<b>Gegenion</b> counter ion		Na <sup>+</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 0.6	0.7
<b>Mindesthaltbarkeit</b> minimum shelf life			11/ 2018
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 13.11.2013

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.