

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>DOWEX® 50 WX2</b> p.A. (100 – 200 mesh), analytical grade	<b>Kat.Nr./Cat.No.</b> : <b>41521</b>
	<b>Lot/Contr.No.</b> : <b>P180202</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance		elfenbeinfarbenes Korn ivory-colored particles	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying	3 h, 105°C	74.0 – 82.0	76.0
<b>Gegenion</b> counter ion		Na <sup>+</sup>	entspricht corresponds
<b>Kapazität (eq/l)</b> capacity		min. 0.6	0.7
<b>Mindesthaltbarkeit</b> minimum shelf life			13.06.2023
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 14.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.