

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

DOWEX® 50 WX8 (100 – 200 mesh) pract.	Kat.Nr./Cat.No. : 41620 Lot/Contr.No. : 150676
---	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		beige-grünes Korn beige-green-colored particles	entspricht corresponds
Trocknungsverlust (%) loss on drying	3 h, 105 °C	50.0 – 58.0	51.4
Gegenion counter-ion		H ⁺	entspricht corresponds
Kapazität (eq/l) capacity		min. 1.7	2.0
Mindesthaltbarkeit minimum shelf life			25.05.2019
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.