

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>SERVA Tris-Glycin/LDS Probenpuffer (4x)</b>	<b>Kat.Nr./Cat.No. : 42525</b>
SERVA Tris-Glycine/LDS sample buffer (4x)	<b>Lot/Contr.No. : P200141</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		rotviolette Lösung red violet solution	entspricht corresponds
<b>pH</b>	1x konzentriert 1x concentrated	7.6 – 7.8	7.8
<b>Zusammensetzung</b> composition		entspricht corresponds	entspricht corresponds
Glycerin/ glycerol		40%	
Lithiumdodecylsulfat lithium dodecyl sulfate		4%	
Ficoll® 400		4%	
EDTA		1.89 mmol/ l	
Bromphenolblau bromophenol blue		0.025%	
Phenolrot/ phenol red		0.025%	
Triethanolamin gepuffert triethanolamine buffered			
<b>Mindesthaltbarkeit</b> minimum shelf life			21.04.2022
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**
**Druckdatum/ printing date: 23.04.2020**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.