

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

<b>SERVA Tris-Glycin/LDS Probenpuffer (4x)</b>	<b>Kat.Nr./Cat.No. : 42525</b>
<b>SERVA Tris-Glycine/LDS sample buffer (4x)</b>	<b>Lot/Contr.No. : P200141</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance		rotviolette Lösung red violet solution	entspricht corresponds
<b>pH</b>	1x konzentriert 1x concentrated	7.6 – 7.8	7.7
<b>Zusammensetzung</b> composition		entspricht corresponds	entspricht corresponds
Glycerin/ glycerol		40%	
Lithiumdodecylsulfat lithium dodecyl sulfate		4%	
Ficoll® 400		4%	
EDTA		1.89 mmol/ l	
Bromphenolblau bromophenol blue		0.025%	
Phenolrot/ phenol red		0.025%	
Triethanolamin gepuffert triethanolamine buffered			
<b>Mindesthaltbarkeit</b> minimum shelf life			21.04.2024
<b>Lagerung (°C)</b> storage			+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 09.05.2022**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.