

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

SERVA Tris-Glycin/LDS Probenpuffer (4x)	Kat.Nr./Cat.No. : 42525
SERVA Tris-Glycine/LDS sample buffer (4x)	Lot/Contr.No. : P230374

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		rotviolette Lösung red violet solution	entspricht corresponds
pH	1x konzentriert 1x concentrated	7.6 – 7.8	7.8
Zusammensetzung composition		entspricht corresponds	entspricht corresponds
Glycerin/ glycerol		40%	
Lithiumdodecylsulfat lithium dodecyl sulfate		4%	
Ficoll® 400		4%	
EDTA		1.89 mmol/l	
Bromphenolblau bromophenol blue		0.025%	
Phenolrot/ phenol red		0.025%	
Triethanolamin gepuffert triethanolamine buffered			
Mindesthaltbarkeit minimum shelf life			16.10.2025
Lagerung (°C) storage			+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.10.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.