

**ANALYSENZERTIFIKAT**  
 CERTIFICATE OF ANALYSIS

<b>Laemmli Probenpuffer (2x), für SDS PAGE</b>	<b>Kat.Nr./Cat.No. : 42526</b>
<b>Laemmli sample buffer (2x), for SDS PAGE</b>	<b>Lot/Contr.No. : P160405</b>

Parameter parameter	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance	blaue Lösung blue solution	entspricht corresponds
<b>pH</b>	6.6 – 7.0	7.0
<b>Zusammensetzung</b> composition		
TRIS/HCl pH 6.8 (# 37180)	126 mM	entspricht corresponds
Glycerin / glycerol (# 23176)	20 %	entspricht corresponds
SDS (# 20763)	4 %	entspricht corresponds
Bromphenolblau (# 15375) bromphenol blue	0.02 %	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		07.11.2018
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 15.11.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.