

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

SERVA Tris-Glycin nativer Probenpuffer (2x)	Kat.Nr./Cat.No. : 42528
SERVA Tris-Glycine native sample buffer (2x)	Lot/Contr.No. : P200270

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	blaue Lösung blue solution	entspricht corresponds
pH	6.6 – 7.0	7.0
Zusammensetzung composition		
TRIS/HCl pH 6.8 (# 37180)	126 mM	entspricht corresponds
Glycerin (# 23176) glycerol	20 %	entspricht corresponds
Bromphenolblau (# 15375) bromphenol blue	0.02 %	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		01.08.2022
Lagerung (°C) storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH Qualitätskontrolle/ Quality Control	Druckdatum/ printing date: 04.08.2020
---	--

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.