

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>SERVA Tris-Glycin/SDS Elektrophorese Puffer (10x)</b>	<b>Kat.Nr./Cat.No. :</b>	<b>42529</b>
<b>SERVA Tris-Glycine/SDS electrophoresis buffer (10x)</b>	<b>Lot/Contr.No. :</b>	<b>P230106</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance	klare, farblose Lösung clear, colorless solution	entspricht corresponds
<b>pH</b>	8.4 – 8.9	8.8
<b>Zusammensetzung</b> composition		
TRIS (# 37180)	0.25 M	entspricht corresponds
Glycin/ Glycine (# 23390)	1.92 M	entspricht corresponds
SDS (# 20765)	1 %	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		03.03.2025
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 13.03.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.