

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

**SERVAGel(TM) IEF Running Buffer Kit**

**Kat.Nr./Cat.No. : 42539**

**Lot/Contr.No. : P150447**

| <b>Parameter</b><br>parameter                   | <b>Spezifikation</b><br>specification                                                         | <b>Ergebnis</b><br>result |
|-------------------------------------------------|-----------------------------------------------------------------------------------------------|---------------------------|
| <b>Anwendung</b><br>application                 | Für die isoelektrische Fokussierung von Proteinen<br>For the isoelectric focusing of proteins | entspricht<br>corresponds |
| <b>Inhalt</b><br>content                        | SERVA IEF Anode Buffer, powder<br>SERVA IEF Cathode Buffer, powder                            | entspricht<br>corresponds |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life |                                                                                               | 30.06.2016                |
| <b>Lagerung (°C)</b><br>storage                 |                                                                                               | +2 bis/ to +8             |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 26.01.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.