

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>SERVAGel(TM) IEF Running Buffer Kit</b>	<b>Kat.Nr./Cat.No. : 42539</b>
	<b>Lot/Contr.No. : P160235</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	Für die isoelektrische Fokussierung von Proteinen For the isoelectric focusing of proteins	entspricht corresponds
<b>Inhalt</b> content	SERVA IEF Anode Buffer, powder SERVA IEF Cathode Buffer, powder	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		07.07.2017
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 10.06.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.