

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS**SERVAGel(TM) IEF Running Buffer Kit****Kat.Nr./Cat.No. : 42539****Lot/Contr.No. : P170300**

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Für die isoelektrische Fokussierung von Proteinen For the isoelectric focusing of proteins	entspricht corresponds
Inhalt content	SERVA IEF Anode Buffer, powder SERVA IEF Cathode Buffer, powder	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		07.07.2018
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 07.09.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.