

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

| | |
|--|--------------------------------|
| TAE Puffer (10x) | Kat.Nr./Cat.No. : 42553 |
| TAE buffer (10x) molecular biology grade | Lot/Contr.No. : P200332 |

| Parameter parameter | Spezifikation specification | Ergebnis result |
|--|---|---------------------------|
| Aussehen appearance | klare, farblose Lösung clear, colorless solution | entspricht corresponds |
| pH | 8.3 – 8.7 | 8.3 |
| Zusammensetzung composition | | |
| TRIS (# 37180) | 48.46 g/l | entspricht corresponds |
| EDTA-Na ₂ ·2H ₂ O (# 11280) | 3.72 g/l | entspricht corresponds |
| Eisessig glacial acetic acid | 12.01 g/l | entspricht corresponds |
| DNase/ RNase Aktivität DNase/ RNase activity | nicht nachweisbar none detected | entspricht corresponds |
| Mindesthaltbarkeit minimum shelf life | | 17.03.2022 |
| Lagerung (°C) storage | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.09.2020

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.