

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

„Laemmli“ Puffer 10x	<b>Kat.Nr./Cat.No. : 42556</b>
„Laemmli“ Buffer 10x	<b>Lot/Contr.No. : P160402</b>

Parameter parameter	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance	klare, farblose Lösung clear, colorless solution	entspricht corresponds
<b>pH</b>	8.4 – 8.9	8.7
<b>Zusammensetzung</b> composition		
TRIS (# 37180)	0.25 M	entspricht corresponds
Glycin / Glycine (# 23390)	1.92 M	entspricht corresponds
SDS (# 20770)	1 %	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		03.11.2018
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 08.11.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.