

**CERTIFICATE OF ANALYSIS**

<b>Semi-dry Blotting Buffer Kit</b> for Western Blotting	<b>Cat.No. :</b> 42559
	<b>Contr.No. :</b> P130517

**Kit Components**

Buffer I	conc. anode buffer	0.3 M Tris and 20 % methanol in aqueous solution
Buffer II	diluted anode buffer	0.03 M Tris and 20 % methanol in aqueous solution
Buffer III	cathode buffer	0.025 M Tris/HCl pH 9.4, 0.04 M 6-aminocaproic acid and 20 % methanol in aqueous solution

**Application**

Western Blotting

**Minimum shelf life**

12/ 2018

**Storage (°C)**

+15 to +30

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date:** 13.12.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.