

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

BlueBlock PF (10x)	Kat.Nr./Cat.No. : 42591
BlueBlock PF (10x)	Lot/Contr.No. : P170141

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		leicht gelbliche Lösung slightly yellowish solution	entspricht corresponds
pH	1:10 verdünnt/ diluted in H ₂ O	7.1 – 7.5	7.3
Anwendung application		entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			06.04.2019
Lagerung (°C) storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 10.04.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.