

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

BlueClear SB, für das Strippen von Antikörpern	Kat.Nr./Cat.No. : 42599
BlueClear SB, for antibody stripping	Lot/Contr.No. : P250122

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Lösung clear, colourless solution	entspricht corresponds
pH		2.6 – 2.8	2.7
Mindesthaltbarkeit minimum shelf life			03.04.2026
Lagerung (°C) storage			+15 bis/ to +30

SERVA Electrophoresis GmbH

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.04.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.