

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Xpress Blotting Puffer (10x) für Western Blotting	Kat.Nr./Cat.No.: 42661
Xpress Blotting Buffer for Western Blotting	Lot/Contr.No.: P180095

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Lösung clear, colourless solution	entspricht corresponds
pH	1:10 verdünnt in H ₂ O 1:10 diluted in H ₂ O	9.1 – 9.5	9.5
Mindesthaltbarkeit minimum shelf life			21.03.2020
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.03.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.