

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

| | |
|--|--------------------------------|
| Xpress Blotting Puffer (10x) für Western Blotting | Kat.Nr./Cat.No. : 42661 |
| Xpress Blotting Buffer for Western Blotting | Lot/Contr.No. : P200084 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|---|--|---------------------------|
| Aussehen appearance | | klare, farblose Lösung clear, colourless solution | entspricht corresponds |
| pH | 1:10 verdünnt in H ₂ O 1:10 diluted in H ₂ O | 9.1 – 9.5 | 9.4 |
| Mindesthaltbarkeit minimum shelf life | | | 27.02.2022 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 02.03.2020

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.