

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Xpress Blotting Kit for Western Blotting	Kat.Nr./Cat.No.: 42662 Lot/Contr.No.: P150417
--	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
-------------------------------	--------------------------	---------------------------------------	---------------------------

Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.
------------------------------------	---

Kitkomponenten
kit-components

250 ml SERVA Xpress Blotting Buffer (10x)

10 Connection Paper (80 x 85 mm)

20 Blotting Fleece sheets (80 x 85 mm)

Mindesthaltbarkeit minimum shelf life	27.11.2016
---	------------

Lagerung (°C) storage	+15 bis/ to +30
---------------------------------	-----------------

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date:** 06.11.2015

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.