

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

---

<b>Xpress Blotting Kit</b> for Western Blotting	<b>Kat.Nr./Cat.No.:</b> 42662 <b>Lot/Contr.No.:</b> P150451
--	--

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
-------------------------------	--------------------------	---------------------------------------	---------------------------

<b>Beschreibung</b> description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.
------------------------------------	---

**Kitkomponenten**  
kit-components

250 ml SERVA Xpress Blotting Buffer (10x)

10 Connection Paper (80 x 85 mm)

20 Blotting Fleece sheets (80 x 85 mm)

<b>Mindesthaltbarkeit</b> minimum shelf life	27.11.2016
---	------------

<b>Lagerung (°C)</b> storage	+15 bis/ to +30
---------------------------------	-----------------

<b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b> We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.
---

---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date:** 01.12.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.