

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS**Xpress Blotting Kit**
for Western Blotting**Kat.Nr./Cat.No.:** 42662
Lot/Contr.No.: P160151

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
-------------------------------	--------------------------	---------------------------------------	---------------------------

Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.		
------------------------------------	---	--	--

Kitkomponenten
kit-components

250 ml SERVA Xpress Blotting Buffer (10x)

10 Connection Paper (80 x 85 mm)

20 Blotting Fleece sheets (80 x 85 mm)

Mindesthaltbarkeit
minimum shelf life

30.08.2017

Lagerung (°C)
storage

+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date:** 05.09.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.