

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS**Xpress PVDF Blotting Kit**  
for Western Blotting**Kat.Nr./Cat.No. : 42664**  
**Lot/Contr.No. : P150446**

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
-------------------------------	--------------------------	---------------------------------------	---------------------------

<b>Beschreibung</b> description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.		
------------------------------------	---	--	--

**Kitkomponenten**  
kit-components

SERVA Xpress Blotting Buffer (10x)

Connection Paper (80 x 85 mm)

Blotting Fleece sheets (80 x 85 mm)

PVDF 0.2 transfer membrane  
(Pore size: 0.2 µm, 80 x 85 mm)**Mindesthaltbarkeit**  
minimum shelf life

27.08.2018

**Lagerung (°C)**  
storage

+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.****Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date: 20.07.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.