

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Xpress PVDF Blotting Kit</b> for Western Blotting	<b>Kat.Nr./Cat.No. : 42664</b> <b>Lot/Contr.No. : P180029</b>
---	--

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
-------------------------------	--------------------------	---------------------------------------	---------------------------

<b>Beschreibung</b> description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.		
------------------------------------	---	--	--

<b>Kitkomponenten</b> kit-components	<p>SERVA Xpress Blotting Buffer (10x)</p> <p>Connection Paper (80 x 85 mm)</p> <p>Blotting Fleece sheets (80 x 85 mm)</p> <p>PVDF 0.2 transfer membrane (Pore size: 0.2 µm, 80 x 85 mm)</p>		
---	---	--	--

<b>Mindesthaltbarkeit</b> minimum shelf life	12.01.2019
---	------------

<b>Lagerung (°C)</b> storage	+15 bis/ to +30
---------------------------------	-----------------

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b> We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 14.02.2018**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.