

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Xpress PVDF Blotting Kit
for Western Blotting

Kat.Nr./Cat.No.: 42664
Lot/Contr.No.: P180216

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
-------------------------------	--------------------------	---------------------------------------	---------------------------

Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.		
------------------------------------	---	--	--

Kitkomponenten
kit-components

SERVA Xpress Blotting Buffer (10x)

Connection Paper (80 x 85 mm)

Blotting Fleece sheets (80 x 85 mm)

PVDF 0.2 transfer membrane
(Pore size: 0.2 µm, 80 x 85 mm)

Mindesthaltbarkeit
minimum shelf life

14.11.2019

Lagerung (°C)
storage

+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 29.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.