

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS



Xpress PVDF Blotting Kit for Western Blotting	Kat.Nr./Cat.No. : 42664 Lot/Contr.No. : P240327
---	--

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
-------------------------------	--------------------------	---------------------------------------	---------------------------

Beschreibung description	Kit für schnelles Semi-Dry Western Blotting von 10 vertikalen mini SDS PAGE Gelen. Kit for fast Semi-Dry Western Blotting of 10 vertical mini SDS PAGE gels.
------------------------------------	---

Kitkomponenten kit-components

SERVA Xpress Blotting Buffer (10x)

Connection Paper (80 x 85 mm)

Blotting Fleece sheets (80 x 85 mm)

PVDF 0.2 transfer membrane
(Pore size: 0.2 µm, 80 x 85 mm)

Mindesthaltbarkeit minimum shelf life

28.03.2026

Lagerung (°C) storage

+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.10.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.