

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

NetFix® for PAG	Kat.Nr./Cat.No. :	42775
265 x 125 mm	Lot/Contr.No. :	P150410

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Größe (mm) size		265 x 125	
Anwendungstest tested for use	Haftungstest adhesion test	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			22.10.2018
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.05.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.