

ANALYSENZERTIFIKAT  
CERTIFICATE OF ANALYSIS



<b>NetFix<sup>®</sup> for PAG</b>	<b>Kat.Nr./Cat.No. : 42775</b>
<b>265 x 125 mm</b>	<b>Lot/Contr.No. : P230080</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Größe (mm)</b> size		265 x 125	
<b>Fadendicke/ Durchmesser</b> Thread thickness/ Yarn diameter		40 µm	entspricht corresponds
<b>Anwendungstest</b> tested for use	Haftungstest adhesion test	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			23.02.2028
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 24.05.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.