

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS



| | | |
|---------------------------------------|----------------------------|----------------|
| SERVALYT® PRECOTES® CSF Kit | Kat. Nr. /Cat.No. : | 42800 |
| for IEF of Cerebrospinal Fluid | Lot/Contr.No. : | P250094 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--|--|---------------------------|
| Beschreibung description | Gebrauchsfertige Polyacrylamidgele mit Servalyt, die an eine Polyesterfolie gebunden sind. Gelformat: 245x125mm, 300µm Ready-to-use polyacrylamide gels with servalyte bound to polyester sheet. Gel format: 245x125mm, 300µm | | |
| pH-Servalyt | 2% in H ₂ O, RT | 6.9 ± 0.3 | 6.90 |
| Leitfähigkeit (mS/cm) conductivity | 2% in H ₂ O, RT | 0.400 – 0.890 | 0.723 |
| Proteintrennung protein separation | IEF, 300 µm-Gel, 5°C, 3500 Vh | entspricht Standard conform reference | entspricht conform |
| | getestet mit CSF-Proben und anschließender Silberfärbung tested with CSF samples and silver staining | | entspricht corresponds |
| Haftung des Gels adhesion of the gel | IEF, 300 µm-Gel, 5°C, 3500 Vh | Haftung ja adhesion yes | ja yes |
| Mindesthaltbarkeit minimum shelf life | SERVALYT® PRECOTES® CSF Gel | | 26.02.2027 |
| Mindesthaltbarkeit minimum shelf life | SERVALYT® PRECOTES® CSF Kit | | 07.05.2026 |
| Lagerung (°C) storage | | | +2 bis/ to +8 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 17.03.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.